

Real-Time PCR Diagnostic Test of Rat Parvovirus Infections and Genetic Strain Identification – Comparative study with serological patterns

INTRODUCTION


The serological screening for Parvovirus infections in rat, has not been quite satisfactory for few years. Mainly, it does not frequently allow to ascertain the virus type, especially for the newly recognized ones. However, this identification may have a major importance, because of their different clinical consequences and interferences with experimental results.

In this study, 38 batches of sentinel rats have been tested, both with serological and molecular biology assays. The comparison of these results is shown hereafter.

MATERIALS AND METHODS


PCR (Scanelis)

Nucleic acids extraction from Rat Caeca Samples (2 to 4 caeca pooled)
Protocol of « High Pure PCR Template Purification Kit », Roche®



Serology (Vebiotel and subcontractors)


- Sampling of caeca and sera was performed by Vebiotel Laboratory for each tested animal
- Serological analyses were performed in OM Diagnostics (The Netherlands) and CRF Laboratories (France)
- Vebiotel analyzed every serological pattern and gave conclusions or hypothesis



SUMMARY

The Rodent Parvovirus group currently includes four different virus types : Kilham Rat Virus (KRV), Toolan's H1 (H1), Rat Minute Virus (RMV) and Rat ParvoVirus (RPV). The virus type involved in an infection must be determined to assess clinical consequences and possible interferences on scientific studies. As there are cross-reactions in serologic assays, it is difficult to type the virus, especially for recently characterized RPV and RMV. However, the association of Real-Time PCR technology and genetic typing by sequencing offers a sensitive and specific diagnostic, especially useful to ascertain seropositive results and to solve inconclusive serologic pattern. Furthermore, acute strain identification allows contaminations monitoring in facilities.

RESULTS 11 positive batches, 2 of them (from 1 institution) infected with KRV and 9 (from 3 institutions) with RMV type



Rodent Parvovirus Real-Time PCR test in Scanelis

- no +/- result and a semi-quantitative approach: sample B has a higher viral load than sample A (about 1000 times more)
- a rodent parvovirus specific and sensitive test


Batches serological positive results compared with molecular typing results in four different institutions

| Batch N° | Sample N° | Facility | KRV ELISA | KRV HAI | KRV IF | H1 ELISA | Serological findings | PCR and sequencing typing (if performed) |
|----------|-----------|----------|-----------|---------|--------|----------|--|--|
| 1 | 1a | B | - | 1:20 | + | - | Parvovirus ? Seroconversion beginning ? | KRV |
| | 1b | | - | - | - | - | KRV | |
| | 1c | | - | - | - | - | KRV | |
| 2 | 2a | A2 | - | - | - | - | RMV ? | RMV |
| | 2b | | - | - | - | - | positive | positive |
| | 2c | | - | - | - | - | positive | positive |
| 3 | 3a | A2 | - | - | - | - | RMV ? | RMV |
| | 3b | | - | - | - | - | positive | positive |
| | 3c | | - | - | - | - | positive | positive |
| 4 | 4a | C | - | 1:80 | + | - | KRV ? RMV ? Both ? | RMV |
| | 4b | | - | - | - | - | RMV | |
| | 4c | | - | - | - | - | RMV | |
| 5 | 5a | A2 | - | - | - | - | RMV ? | RMV |
| | 5b | | - | - | - | - | RMV | |
| | 5c | | - | - | - | - | RMV | |
| 6 | 6a | A2 | - | 1:80 | + | - | KRV ? RMV ? Both ? | RMV |
| | 6b | | - | - | - | - | RMV | |
| | 6c | | - | - | - | - | RMV | |
| 7 | 7a | A2 | - | 1:80 | + | - | RMV ? | positive |
| | 7b | | - | - | - | - | RMV ? | positive |
| | 7c | | - | - | - | - | RMV ? | positive |
| 8 | 8a | A1 | - | - | - | - | RMV ? | negative |
| | 8b | | - | - | - | - | RMV ? | RMV |
| | 8c | | - | - | - | - | RMV ? | positive |
| 9 | 9a | A1 | - | - | - | - | RMV ? | RMV |
| | 9b | | - | - | - | - | RMV ? | positive |
| | 9c | | - | - | - | - | RMV ? | positive |
| 10 | 10a | D | - | 1:20 | + | - | RMV ? | positive |
| | 10b | | - | - | - | - | RMV ? | positive |
| | 10c | | - | - | - | - | RMV ? | positive |
| 11 | 11a | D | - | - | - | - | Seroconversion beginning ? | negative |
| | 11b | | - | - | - | - | Seroconversion beginning ? | negative |

| | Seropositive | Seronegative |
|--------------|--------------|--------------|
| PCR positive | 11 | 0 |
| PCR negative | 1 +/- sero | 26 |


Serological versus PCR result for each batch

- PCR and serological results provide a good agreement for one batch
- Whereas the serological patterns inconsistently allow viral identification, the association of PCR and sequencing provides the Rat Parvovirus type.
- According to these results, a limited number of sentinels may be tested with PCR assay
- PCR enables to confirm or infirm the infection, when serological patterns are inconclusive



Comparative Analysis of RMV isolates DNA sequences in one institution:

Genetic stability of a Rat Minute Virus isolate during a relative long period in a same facility and Rat Minute Virus typing in a same institution for contamination monitoring



Legend:
 ◊ ◊ ◊ different nucleotides (A, T, G et C)
 Total homology (100%)
 Facility 1
 Facility 2
 Colours in symbols mean one base mutation between the different RMV isolates

- 100% homology over 4 months seems to mean a relative genetic stability of this RMV isolate
- Molecular typing: a tool for the approach to the RMV's strain. The variability of isolates in a same institution helps to implement for contamination management
- The 162 pb product sequencing is sufficient to determine the Rat Parvovirus Type. However, it would be interesting to sequence a longer fragment of the viral genome, to obtain a more accurate identification of the viral isolate. It may improve the contamination monitoring in an institution

CONCLUSIONS

- Serological and PCR assays are complementary for an early and precise diagnosis of parvoviral infections
- This molecular assay allows :
 - an accurate diagnosis and typing of the parvovirus species : KRV, Toolan H1, RPV (RPV-1a) and RMV (RPV2)
 - confirmation or infirmation of an inconclusive serological pattern

Perspectives:

- Environmental monitoring
- Screening in non-immunocompetent rats
- Testing others samples may enhance sensitivity (mesenteric lymph node) or allow diagnosis on live animal (faeces) ?
- Pooling samples may reduce analysis cost